



System products, nomenclature

- All system products are supplied in an individual sterile pouch.
- Avance® Solo and Avance® Solo Adapt are supplied in a modular format to allow the user to select according to their requirements and to minimise potential for overstocking or waste.
- Individual sterile units; products combined with a single sterile pouch or blister pack.
- Retail pack: a retail pack may contain more than one individual sterile unit.

Avance® Solo Assortment*

Product Code	Description / Size (cm)	Individual Sterile Unit Contents	Pcs per box
Pump, canister and border dressing			
8881020	Avance Solo Pump, Canister and Border Dressing - 10 x 20cm	Pump + Canister 50ml with tubing + 4 x AA batteries + 10 x 20cm dressing + Belt clip	1
8881030	Avance Solo Pump, Canister and Border Dressing - 10 x 30cm	Pump + Canister 50ml with tubing + 4 x AA batteries = 10 x 30cm dressing + Belt clip	1
8881035	Avance Solo Pump, Canister and Border Dressing - 10 x 35cm	Pump + Canister 50ml with tubing + 4 x AA batteries + 10 x 35cm dressing + Belt clip	1
Avance Solo pump and canister			
8880050	Avance Solo Pump and Canister - 50ml	Pump + Canister 50ml with tubing + 4 x AA batteries + Belt clip	1
Border dressing			
881020	Avance Solo Border Dressing - 10 x 20cm	10 x 20cm Border dressing + 6 secondary fixation strips	2
881030	Avance Solo Border Dressing - 10 x 30cm	10 x 30cm Border dressing + 6 secondary fixation strips	2
881035	Avance Solo Border Dressing - 10 x 35cm	10 x 35cm Border dressing + 6 secondary fixation strips	2
881515	Avance Solo Border Dressing - 15 x 15cm	15 x 15cm Border dressing + 6 secondary fixation strips	2
881520	Avance Solo Border Dressing - 15 x 20cm	15 x 20cm Border dressing + 6 secondary fixation strips	2
881530	Avance Solo Border Dressing - 15 x 30cm	15 x 30cm Border dressing + 6 secondary fixation strips	2
882020	Avance Solo Border Dressing - 20 x 20cm	20 x 20cm Border dressing + 6 secondary fixation strips	2
882525	Avance Solo Border Dressing - 25 x 25cm	25 x 25cm Border dressing + 6 secondary fixation strips	2
Other system consumables			
880050	Avance Solo Canister 50ml	Canister 50ml with tubing	4
882000	Avance Solo Foam 1.5 x 10 x 12cm	Foam 1.5 x 10 x 12cm	4
Other system consumables			
890001	Carry bag small	15 x 12 x 3.5cm	7
890002	Carry bag large	30.5 x 14.5 x 7.5cm	7
5105100-00	Nightcase		8

Avance® Solo Adapt Assortment*

Product Code	Description / Size (cm)	Individual Sterile Unit Contents	Pcs per box
9994152	Avance Solo Adapt Pump, Canister, Foam, Film and Transfer Port	Pump + Canister 50ml with tubing + 4 x AA batteries + Safetac film dressing + Foam (8 x 10 x 3cm) + Belt clip	1
994152	Avance Solo Adapt Foam, Film and Transfer Port	Safetac film dressing + Foam (8 x 10 x 3cm) + Transfer port and tubing	5
994000	Avance Solo Adapt Film - 20 x 27cm	Safetac adhesive film dressing	5

www.molnlycke.com.au | www.molnlycke.co.nz

You will be required to input the following code **1FYrV** or your **AHPRA registration number** to access our websites.

Mölnlycke, 12 Narabang Way, Level 4, Betrose, NSW, 2085. Phone 1800 005 231.
New Zealand Orders & Enquiries 0800 005 231. www.molnlycke.co.nz. The Mölnlycke, Avance and Safetac trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2022 Mölnlycke Health Care AB. All rights reserved. AUWC_000998



Avance® Solo and Avance® Solo Adapt

Product descriptions and units of supply



Avance® Solo

– Closed incisions and open wounds (ciNPT and dNPWT)


1 Pump, canister and border dressing

The three sets available in the Avance Solo assortment are designed for convenience in the OR primarily for closed incision applications. These sets contain everything required for an application and commencement of the therapy without duplications so as to avoid waste. A pump and a canister plus a dressing, fixation strips in one package and 4 x AA Li batteries.



3 Border dressing

The eight different sizes of dressings are supplied in retail packaging containing two individually packaged sterile dressings. Each dressing has an affixed dressing tube and a male connector, six polyurethane film secondary fixation strips with an acrylic adhesive.



2 Avance Solo pump and canister

The Avance Solo pump is always packaged in a format that includes:

1. Pump body housing the motor and pump assembly
2. A plastic belt clip
3. 50ml canister with tubing and a female connector
4. 4 x AA Li batteries



4 Avance Solo canister and foam

Canisters and a foam wound filler are supplied in individually pouched and sterilised pieces that are combined as a group of four in a single package. Canisters have pre-affixed tubing and a connector.




Avance® Solo Adapt

– Pressure injuries

1 Avance Solo Adapt

The Avance® Solo Adapt system uses the same canister as Avance Solo, but the pump is coupled to a Safetac adhesive film dressing that can be cut and customised to support an application over difficult anatomical areas, such as the sacral region.



2 Avance Solo Adapt Dressing

Dressings are supplied in an individual sterile pouch containing a Safetac® adhesive film dressing, Avance Solo Adapt Foam, tubing and a connector.



3 Safetac Film Dressing

Safetac adhesive film dressings are supplied in an individual sterile pouch.

